

# Super Deluxe

Hydraulic General OT Table



# POSITIONS

Chair Position



Lithotomy Position



Trendelenburg/  
Anti-Trendelenburg Position



Lateral Tilt Position



Height Adjustment  
Position



Kidney Position



## Technical Constraints

<b>Top dimension</b>	<b>L 1854 x W 500mm</b>
<b>Height adjustment</b>	<b>750mm- 1000mm</b>
<b>Trendelenburg/ Reverse</b>	<b>30° / 25°</b>
<b>Lateral tilt</b>	<b>20° / 20°</b>
<b>Kidney elevator</b>	<b>150mm</b>
<b>Back Rest (up/down)</b>	<b>80°/ 25°</b>
<b>Leg Rest (up/down)</b>	<b>15° / 90°</b>
<b>Head Rest (up/down)</b>	<b>20° / 60°</b>
<b>Patient Weight Capacity</b>	<b>120kgs (264lbs)</b>

## STANDARD ACCESSORIES



**Knee/Lithomic Crutches**



**Lateral Support**



**Anaesthesia Screen**



**Arm Board**



**Shoulder Support**

---

**For Further Information, contact:**

**# 36 B, Bajrang Bali Colony, Rohtak Gate,  
Bhiwani - 127021 (HR)**

**M.: 9802324555, 8826707666**

**Email: [nvsurgicals666@gmail.com](mailto:nvsurgicals666@gmail.com)**