

# C-Arm Compatible Hydraulic OT Table



# POSITIONS



## Technical Constraints

<b>Top dimension</b>	<b>L 1980 x W 533mm</b>
<b>Height adjustment</b>	<b>750mm-1000mm</b>
<b>Trendelenburg/ Reverse</b>	<b>30° / 25°</b>
<b>Lateral tilt</b>	<b>20° / 20°</b>
<b>Kidney elevator</b>	<b>150mm</b>
<b>Back Rest (up/down)</b>	<b>80° / 25°</b>
<b>Leg Rest (up/down)</b>	<b>15° / 90°</b>
<b>Head Rest (up/down)</b>	<b>20° / 60°</b>
<b>Patient Weight Capacity</b>	<b>150KGS (330LBS)</b>

## STANDARD ACCESSORIES



**Knee/Lithomic Crutches**



**Lateral Support**



**Anaesthesia Screen**



**Arm Board**



**Shoulder Support**

---

**For Further Information, contact:**

**# 36 B, Bajrang Bali Colony, Rohtak Gate,  
Bhiwani - 127021 (HR)**

**M.: 9802324555, 8826707666**

**Email: [nvsurgicals666@gmail.com](mailto:nvsurgicals666@gmail.com)**