

Electromatic C-Arm

Compatible OT Table



POSITIONS



POSITIONS



Technical Constraints

Top dimension	L 1980 x W 533mm
Height adjustment	750mm-1000mm
Trendelenburg/ Reverse	30° / 25°
Lateral tilt	20° / 20°
Kidney elevator	150mm
Back Rest (up/down)	80° / 25°
Leg Rest (up/down)	15° / 90°
Head Rest (up/down)	20° / 60°
Power supply	24VDC
Patient Weight Capacity	200KGS (440LBS)

STANDARD ACCESSORIES



Knee/Lithomic Crutches



Lateral Support



Anaesthesia Screen



Arm Board



Shoulder Support

For Further Information, contact:

**# 36 B, Bajrang Bali Colony, Rohtak Gate,
Bhiwani - 127021 (HR)**

M.: 9802324555, 8826707666

Email: nvsurgicals666@gmail.com